

# #9200 CHILDREN'S VECTOGRAPHIC PROJECTOR SLIDE INSTRUCTIONS

The Stereo Optical Children's Vectographic Test Slide has been designed for use with Marco, A.O. (Reichert), Woodlyn and similar projectors.

## CONVERT FOR USE IN A 10 FOOT LANE

Although initially guredcon for use in a 20 foot lane, the 9200 test slide can be easily converted for a 10 foot lane.

1. The glass slide is held in the frame by a spring in the bottom channel.
2. Remove the glass slide from the frame as follows:
  - ◆Grasp slide with one hand and frame with the other.
  - ◆Slide glass down within frame.
  - ◆Tilt top of slide out and remove.
3. Rotate glass slide on a vertical axis.
4. Return slide to frame by reversing process step 2.

## ADMINISTERING THE VECTOGRAPHIC TESTS

The 9200 slide will administer both standard and vectographic testing - there is no need to change the slide.

The following tests must be viewed under polarized conditions:

- ◆Binocular Balance
- ◆Suppression
- ◆Fixation Disparity
- ◆Stereopsis

These tests evaluate the performance of each eye individually, while both eyes remain open and fused.

Polarized lters can be purchased through Stereo Optical or one of its distributors. For the distributor nearest you, call 1-800-344-9500.

## CLEANING INSTRUCTIONS

**Do not** immerse slide in any liquid or spray liquid directly on slide. Use a soft damp cloth to clean.

# #9200 CHILDREN'S VECTOGRAPHIC PROJECTOR SLIDE ANSWER KEY

## LETTER ACUITY (Left to Right)

LEVEL	ANSWER
20/200	Right
20/160	Top Right
20/125	Right Left
20/100	Left Bottom Right
20/80	B R T
20/60	R L T R
20/50	B L T B L
20/40	R B R T L
20/30	T R B --- L B T
20/25	L B B R --- L T R B
20/20	T R L B --- R B T L
20/15	R T L R --- T B L T

## BINOCULAR BALANCE

Ring, Hand, Bat, Plus Sign, Infant

LEVEL	LEFT	RIGHT
20/30	T L B R	R B T L
20/25	B L T B	T L B R
20/20	R T L B	B L R T

## SUPPRESSION

Sequencing: Binocular, Left, Right, Binocular, Left, Right...


20/40	L R R T B L
20/30	B T B L R L T
20/25	R T B L L T R B R
20/20	T B L T R B L R B T
20/15	R L B B R T L T B R

## FIXATION DISPARITY

Left: Arrow Up, Line Left  
Right: Arrow Down, Line Right

## STEREOPSIS (which of the rings stands out?)

LEVEL	ANSWER
Zero	S qua re
1 min of Arc	Circle
2 min of Arc	Cross
3 min of Arc	Triangle
4 min of Arc	Star



## STEREO OPTICAL CO., INC.

8600 W. Catalpa Avenue, Suite 703  
Chicago, IL 60656 - USA  
Phone: 1.800.344.9500  
Email: sales@stereooptical.com  
Website : www.StereoOptical.com



©2017 Stereo Optical Co., Inc. All rights reserved. Unless indicated otherwise, all trademarks are the property of Stereo Optical Co., Inc. in the United States and other countries.