

STEREO OPTICAL DRIVER REHAB RECORD FORM FOR THE OPTEC[®] VISION TESTER

NAME: _____ DATE: _____
 ADDRESS: _____ AGE: _____
 ADDRESS: _____ WEARING GLASSES: Y N
 ADMINISTRATOR: _____ CONTACT LENSES: Y N

PERIPHERAL TEST $L = 85^\circ 70^\circ 55^\circ$ 40° NASAL 40° $R = 85^\circ 70^\circ 55^\circ$

DIAL AT 01: ACUITY - FAR Score: 20/_____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	Z N	R O	H K
2 20/100	R K S	H N C	Z O D
3 20/70	H C D V	S K Z O	R N D S
4 20/50	Z R O D	N S C H	V Z K N
5 20/40	K H S C	O Z N R	D N V C
6 20/30	O N R Z V	D K H C S	K D S O N
7 20/20	S D C H N	V R Z K O	H S N R D

DIAL AT 03: ACUITY - NEAR Score: 20/_____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	S V C	N R K	H Z O
2 20/70	R N Z H	D O K V	C S Z N
3 20/50	C K V D	S N Z R	D O H C
4 20/40	V H R N	O D S K	N Z C S
5 20/30	H S K R C	N Z D O V	Z S H N K
6 20/20	Z O N V R	H C S K D	V K C D S

DIAL AT 02: COLOR PERCEPTION - FAR A-1 2 B-5 C-2 6 D-6 E-1 6 F-BLANK PASS ____ FAIL ____

DIAL AT 04: LATERAL PHORIA - FAR 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15

1 PRISM DIOPTERS E S O P H O R I A O R T H O E X O P H O R I A

DIAL AT 05: VERTICAL PHORIA - FAR 1 2 3 4 5 6 7

1/2 PRISM DIOPTERS L E F T H Y P E R P H O R I A O R T H O R I G H T H Y P E R P H O R I A

DIAL AT 06: LATERAL PHORIA - NEAR 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15

1 PRISM DIOPTERS E S O P H O R I A O R T H O E X O P H O R I A

DIAL AT 07: DEPTH PERCEPTION - FAR 1-B 2-L 3-B 4-T 5-T 6-L 7-R 8-L 9-R

SECONDS OF ARC 400 200 100 70 50 40 30 25 20

DIAL AT 08: FUSION - FAR 3 cubes BOTH eyes = PASS Other than 3 cubes = FAIL

DIAL AT 09: LEFT EYE Score: 20/_____ Right Occluder OFF **RIGHT EYE** Score: 20/_____ Left Occluder OFF

TUMBLING "E" FAR	PILOT LINE	LEFT EYE								RIGHT EYE							
		20		20		20		20		20		20		20		20	
		200	100	70	50	40	30	20	200	100	70	50	40	30	20		
LINE/ACUITY	COLOR	L	T	B	L	R	L	T	LINE/ACUITY	R	B	R	L	L	T	B	
A 20/50	BLUE		R	T	R	L	T	B	A 20/50		L	R	T	L	B	T	
B 20/40	RED		B	R	B	L	L	R	B 20/40		T	T	L	B	T	R	
C 20/30	AMBER		L	T	T	L	R	B	C 20/30		L	B	R	R	L	T	
D 20/20	GREEN		R	L	T	L	B	T	D 20/20		R	T	R	L	T	B	

DIAL AT 10: ACUITY - FAR Score: 20/_____

LINE/ACUITY	COLOR	LEFT	BOTH	RIGHT
1 20/200	RED	3 9	6 5	2 8
2 20/100	GREEN	8 5 2	3 9 5	6 2 8
3 20/70	RED	5 3 6 2	9 8 5 6	2 3 6 5
4 20/50	AMBER	9 2 8 6	5 9 6 3	2 8 5 9
5 20/40	GREEN	5 8 6 9	2 6 8 3	9 5 3 2
6 20/30	RED	3 2 5 9 6	6 8 3 5 2	8 2 6 3 5
7 20/20	AMBER	2 9 8 6 5	3 5 2 6 9	6 5 3 8 2

- DIAL AT 11: SIGN RECOGNITION AND DEPTH PERCEPTION**
1. HILL
 2. HANDICAPPED PARKING
 3. * DIVIDED HIGHWAY (581 Sec of Arc)
 4. NO PARKING
 5. REGULATION
 6. * YIELD (323 Sec of Arc)
 7. NO U TURN
 8. DO NOT ENTER
 9. STOP SIGN
 10. RR CROSSING
 11. * SCHOOL ZONE (145 Sec of Arc)
 12. NO PASSING

DIAL AT 12: COLOR RECOGNITION

LINE A	1=AMBER	2=GREEN	3=RED	4=BLUE
LINE B	1=RED	2=AMBER	3=BLUE	4=GREEN
LINE C	1=GREEN	2=BLUE	3=AMBER	4=RED
LINE D	1=BLUE	2=RED	3=GREEN	4=AMBER

© 2017 Stereo Optical Co., Inc. All rights reserved. Unless indicated otherwise, all trademarks are the property of Stereo Optical Co., Inc. in the United States and other countries.