


# Optec® 1000 Test slide options


Multiple Acuity Slide (1000-23) Standard  
Acuity slide for preschool children and adults



Test 1 & 2 Far Point  
Acuity Allen Test  
Binocular (20/100)  
Monocular (20/50 to 20/30)

Test 3 & 4 Far Point  
Acuity  
Right Eye & Left Eye Monocular  
or Both Eyes Simultaneously  
(20/200 to 20/20)

Elementary Test Slide (1000-22) Standard  
Vision screening slide for elementary school age children




Test 1 & 2 Far Point  
Acuity Right Eye & Left  
Eye Monocular (20/200 to 20/20)  
Also Right Eye & Left Eye Plus Lens Test

Test 3 Far Point  
Muscle Balance  
(Combination of Lateral and Vertical Phoria)

Test 4 Far Point  
Color Perception Test, (Tumbling "E")

Secondary Test Slide (1000-21) Optional  
Acuity and Phoria slide for children of secondary school age and adults



Test 1 & 2 Far Point  
Acuity Right Eye & Left Eye  
Monocular or Both Eyes Simultaneously  
(20/200 to 20/20)

Test 3 Far Point  
Lateral Phoria Test  
(1 Diopter Increments)

Test 4 Far Point  
Vertical Phoria Test  
(1/2 Diopter Increments)